

LimeLight provides the speed and versatility of an IPL with the accuracy and precision of a laser. Cutera's real-time calibration feature, programmable wavelengths, integrated dynamic cooling and consistent pulse maximize both safety and clinical outcomes. In the world of IPLs, LimeLight holds an unique position for clinical excellence.

Practice benefit

- High efficacy and low risk
- No downtime
- Designed with and endorsed by Asian physicians
- Programmable wavelengths
- No disposables



Lentigines

Skin type IV, Japanese, Female

Treatment Protocol

Program: A mode
Cooling: Normal mode (5°C)
Fluence: 12 J/cm²



Before



After 1 Tx

Photos courtesy of Kei Negishi, M.D. Tokyo Women's Medical University Aoyama Institute of Women's Medicine, Tokyo, Japan