

LimeLight provides the speed and versatility of an IPL with the accuracy and precision of a laser. Cutera's real-time calibration feature, programmable wavelengths, integrated dynamic cooling and consistent pulse maximize both safety and clinical outcomes. In the world of IPLs, LimeLight holds an unique position for clinical excellence.

### Practice benefit

- High efficacy and low risk
- No downtime
- Designed with and endorsed by Asian physicians
- Programmable wavelengths
- No disposables



## Lentigines

Skin type IV, Japanese, Female

## Treatment Protocol

Program: B mode  
Cooling: Normal mode (5°C)  
Fluence: 12 J/cm<sup>2</sup>



Before



After 1 Tx

Photos courtesy of Kei Negishi, M.D. Tokyo Women's Medical University Aoyama Institute of Women's Medicine, Tokyo, Japan