

# Dornier *Medilas D LiteBeam+*

1470 nm  
wavelength

 **LPS**



# Dornier *Medilas D LiteBeam+*

With an increasing number of vein disorders among the general population, the ability to treat effectively and efficiently with minimal side effects has become a requirement. Dornier's proven diode technology is the solution for all of your endovenous laser needs. Incorporating the 1470 nm wavelength and an output power of up to 15 Watts, procedures are performed with less patient discomfort using the **Medilas D LiteBeam+**. Unique safety features such as the Lightguide Protection System (LPS) and fibertom mode minimizes risk to the patient and fiber.



Before



After

## **Benefits of LiteBeam+**

- *Less Patient Discomfort*
- *Compact & Transportable*
- *Variety of Flexible Fibers*
- *Intuitive Graphical Interface*
- *Enhanced Safety Features*

## **Treatment Flexibility**

The **Medilas D LiteBeam+** can be used with various sized thin flexible fibers. Extending treatment flexibility, Dornier offers both single and multiple use fibers.

## **fibertom Mode**


This feature produces even cutting by ensuring the temperature of the fiber tip remains within set limits. This provides well controlled application of energy while reducing carbonization and preventing damage to the fiber tip.

## **Ergonomic by Design**

The Dornier **Medilas D LiteBeam+** is lightweight and easily transported using the integrated handles. The compact design of this tabletop laser allows it to be utilized in a variety of clinical and office settings, while the intuitive color touch-screen is easy to use.



## 1470 nm wavelength



### **Gold Standard Treatment**

Endovenous laser treatment has replaced ligation and stripping as the Gold Standard of Care. The benefits of this revolutionary procedure include:

- *Office-based procedure*
- *Local anesthesia*
- *30 – 45 minute procedure*
- *Virtually no post-operative patient recovery time*
- *Minimal learning curve for the physician*
- *Cost effective*
- *Reimbursement available*
- *Increased physician revenue streams*
- *Equipment has FDA-clearance*

### **Lightguide Protection System (LPS)**

Based upon proprietary technology, the Dornier *Medilas D LiteBeam+* comes equipped with an unique Lightguide Protection System (LPS) that protects the fiber from damage thus increasing patient safety. If light emissions originating from excess heat during the application of laser energy is detected at the applicator tip, the LPS system turns off the laser to prevent damage to the fiber.





WE'RE ALL ABOUT PEOPLE.

## Headquarters

Dornier MedTech GmbH  
Argelsrieder Feld 7  
D-82234 Wessling  
Germany

## America

Dornier MedTech America, Inc.  
1155 Roberts Boulevard  
Suite 100  
Kennesaw, GA 30144  
USA  
Phone: +1-770 426-1315  
Fax: +1-770 426 6115  
Email: [info@dornier.com](mailto:info@dornier.com)

## Asia/Pacific

Dornier MedTech Asia Pte Ltd  
2 Corporation Road  
#05-11/12 Corporation Place  
Singapore 618494  
Phone: +65-6665-5788  
Fax: +65-6665 6711  
Email: [infoasia@dornier.com](mailto:infoasia@dornier.com)

## Japan

Dornier MedTech Japan  
Meguro Estate Building  
3-8-5 Kami-Ohsaki  
Shinagawa-ku, Tokio 141-0021  
Japan  
Phone: +81-3-3280 3550  
Fax: +81-3-3280 3555  
Email: [info@domedtech.co.jp](mailto:info@domedtech.co.jp)

## Europe

Dornier MedTech Europe GmbH  
Argelsrieder Feld 7  
D-82234 Wessling  
Germany  
Phone: +49-8153-888 625  
Fax: +49-8153-888 444  
Email: [infoeurope@dornier.com](mailto:infoeurope@dornier.com)

Dornier MedTech Italia s.r.l.  
Via Arrigo Cavaglieri N. 26  
00173 Rome  
Italy  
Phone: +39-06-72 35 289  
Fax: +39-06-72 35 302


Dornier MedTech Espana, S. L.  
Pedro Y Pons 9-11  
Edificio Masters  
E-08034 Barcelona  
Spain  
Phone: +34-93-203 9316  
Fax: +34-93-205 3814

Dornier MedTech France SARL  
Ferme de St Paul  
1545 route d'Epagny  
F - 74330 Sillingy  
France  
Phone: +33-4-50 22 18 94  
Fax: +33-4-50 24 24 53

Dornier MedTech Europe GmbH  
Moscow Representative Office  
Uliza Mytnaya 3  
Office 21  
119049 Moscow  
Russia  
Phone: +7-495-230 29 94  
Fax: +7-495-739 51 33



[www.dornier.com](http://www.dornier.com)

 1275

[info@dornier.com](mailto:info@dornier.com)



DMTA345-A

**Dornier** *MedTech*